


# Importing a Respondus Test Into TurningPoint®

TurningPoint 4.x  
Respondus 3.5

The following FAQs have been provided for your convenience. If you have additional support questions, please feel free to contact Turning Technologies Tech Support at 1-866-746-3015 or visit [www.TurningTechnologies.com](http://www.TurningTechnologies.com).

## To Save a TurningPoint Compatible Test:

1. Open your Respondus software.
2. In the **Start | Open** options, select a **Current Personality** of "IMS QTI 1.1+".
3. Create or Edit your test using Multiple Choice, Multiple Response, True/False, or Fill-in-the-Blank questions.
4. Click the **Preview & Publish** tab.
5. Choose "**Save QTI XML File.**"
6. In the **File Format** option, select "**QTI XML zip file**", then select "**Save As**". Name your file and save it for later import to TurningPoint.



## To Import a Compatible Respondus Test into TurningPoint:

1. Select **Tools** from the TurningPoint toolbar.
2. Select **Parser**.
3. Select **QTI Document**.
4. Select your saved Respondus zip file.
5. TurningPoint will create slides from each of your questions in your Respondus Test.

